

REF

DK.045.01.9



For *In vitro* Diagnostic Use

INSTRUCTIONS FOR USE

abia AFP

**Enzyme immunoassay for the quantitative determination
of alpha-fetoprotein (AFP) concentration
in human serum**

This Package Insert provides information for Professional Use of the kit.

The kit contains sufficient reagents for 96 (breakable wells) assays including controls; the kit is intended for manual testing with a possibility of fractional (one well) use of the kit or use of the kit on open type automated analyzer for enzyme immunoassay.

I. INTENDED USE

The abia AFP kit is intended for the quantitative determination of Alpha-Fetoprotein (AFP) concentration in human serum by a microplate immunoenzymometric assay.

This kit is for diagnostic use by a trained laboratory professional and will not be sold to the general public. All the reagents are for professional *in vitro* diagnostic use only.

The results of this or any other diagnostic assay should be used and interpreted only in the context of the overall clinical picture.

II. INTRODUCTION

Alpha-fetoprotein (AFP) is a glycoprotein with a molecular weight of approximately 70000 Daltons. AFP is normally produced during fetal and neonatal development by the liver, yolk-sac and in small concentrations by the gastrointestinal tract. After birth, serum AFP concentrations decrease rapidly, and by the second year of life and thereafter only trace amounts are normally detected in serum (less than 10 IU/ml).

Elevation of serum AFP to abnormally high values occurs in several malignant diseases, most notably nonseminomatous testicular cancer and primary hepatocellular carcinoma. In the case of nonseminomatous testicular cancer, a direct relationship has been observed between the incidence of elevated AFP levels and the stage of disease. Elevated AFP levels have also been observed in patients diagnosed with seminoma with nonseminomatous elements, but not in patients with pure seminoma.

In addition, elevated serum AFP concentrations have been measured in patients with other noncancerous diseases, including ataxia telangiectasia, hereditary tyrosinemia, neonatal hyperbilirubinemia, acute viral hepatitis, chronic active hepatitis and cirrhosis. Elevated serum AFP concentrations are also observed in pregnant women. Therefore, AFP measurements are not recommended for use as a screening procedure to detect the presence of cancer in the general population.

III. PRINCIPLE OF THE TEST

The abia AFP is a one-step immunoassay, based on principle of “sandwich” method. The assay system utilizes a high affinity and specificity monoclonal antibody (enzyme conjugated and immobilized) directed against a distinct antigenic determinant on the intact AFP molecule. The test sample is allowed to react simultaneously with the two antibodies, resulting in the AFP molecules being sandwiched between the solid phase and enzyme-linked antibodies. After incubation, the wells are washed with Washing Solution to remove unbound labeled antibodies. A solution of TMB-Substrate is added and incubated, resulting in the development of a blue color. The color development is stopped with the addition of Stopping Reagent, changing the color to yellow. The color intensity of the test sample is directly proportional to the concentration of AFP. Absorbance is measured spectrophotometrically at 450 nm.

IV. CONTENT OF THE KIT abia AFP

Table 1

LABEL	NATURE OF THE REAGENTS	PRESENTATION
Anti-AFP-coated microtiter wells	Polystyrene stripped 96-well plate (breakable wells) coated with monoclonal antibodies to AFP. Store at 2-8 °C until expiration date.	1 plate
Conjugate	Monoclonal anti-AFP antibodies conjugated to horseradish peroxidase. Contains bovine serum albumin (5.0%) and Tween-20 (0.1%). Transparent or opalescent pink color liquid. Preserving agent: 0.05% Proclin 300, 0.004% gentamycin sulfate, 0.1% phenol. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 12.0 ml
Calibrator 0 Calibrator 1 Calibrator 2 Calibrator 3 Calibrator 4 Calibrator 5	Six vials containing AFP and bovine serum albumin (1.0%). Calibrators were calibrated using a WHO 1 st IRP 72/225. The AFP concentration levels in Calibrators are provided on the labels of vials and in the Certificate of Analysis on a lot-specific basis. Transparent or opalescent liquids, pale yellow. Preserving agent: 0.05% Proclin 300, 0.004% gentamycin sulfate, 0.1% phenol. Store at 2-8 °C until expiration date in a tightly sealed vial.	5 vials 0.5 ml. Calibrator 0 – 2.0 ml
Control Serum	Human serum with a defined quantity of AFP, containing bovine serum albumin (1.0%). The AFP concentration level in Control Serum is provided on the vial label and in the Certificate of Analysis on a lot-specific basis. Transparent or opalescent liquid, pale yellow. Preserving agent: 0.05% Proclin 300, 0.004% gentamycin sulfate, 0.1% phenol. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 0.5 ml
Washing Solution (concentrated 25-fold)	Phosphate-saline solution (pH 7.4-7.7). Transparent or slightly opalescent liquid, colorless, or pale yellow, sediment may form that dissolves completely at 35-39 °C and shaking. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 50.0 ml
TMB-Substrate	Tetramethylbenzidine (0.03%) in citric acid buffer, containing H ₂ O ₂ (0.01%). Transparent colorless liquid. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 14.0 ml
Stopping Reagent	0.2M sulfuric acid solution. Transparent colorless liquid. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 25.0 ml
Protective films for EIA plates		2
Disposable tips		16
Disposable plastic dishes for liquid reagents		2
Polyethylene bag with a Zip-Lock		1

V. PRECAUTIONS

The reliability of the results depends on correct implementation of the following requirements:

- The temperature in the lab should be 18-24 °C.
- Inspect the contents of the box: check the vials and labels integrity. In case of label loss or labels/vials damage, vials should be disposed and kit cannot be used.
- Do not use expired reagents.
- Do not mix reagents from different lots within a given test run.
- Carefully reconstitute the reagents avoiding any contamination.
- Do not carry out the test in the presence of reactive vapors (acid, alkaline, aldehyde vapors) or dust that could alter the enzyme activity of the conjugates.
- Use glassware thoroughly washed and rinsed with deionized water or preferably, disposable material.

- Do not allow the microplate to dry between the end of the washing operation and the reagent distribution.
- The enzyme reaction is very sensitive to metal ions. Consequently, do not allow any metal element to come into contact with the various conjugate or TMB-Substrate.
- Use a new distribution tip for each sample.
- Do not reuse protective films for EIA plates.
- Well washing is a critical step in this procedure: respect the recommended number of washing cycles and make sure that all wells are completely filled and then completely emptied. Incorrect washing may lead to inaccurate results.
- Never use the same container to distribute conjugate and other solutions.
- Check the pipettes and other equipment for accuracy and correct operation.
- Do not change the assay's procedure.
- Use distilled or deionized water.
- Avoid exposure of the reagents to excessive heat or sunlight during storage and incubation.
- Once the assay has been started, all subsequent steps should be performed without interruption. Do not let the wells dry once the assay has been started.

VI. HEALTH AND SAFETY INSTRUCTIONS

- All reagents included in the kit are intended for “in vitro diagnostic use”.
- Human origin material used in the preparation of Control Serum and Calibrators have been tested and found negative for HBsAg, antibodies to hepatitis C virus and antibodies to human immunodeficiency virus (HIV-1 and HIV-2).
- Certain reagents contain biological material of animal origin.
- Because no known test method can offer complete assurance that infectious agents are absent, handle reagents and patients samples as if capable of transmitting infectious disease.
- Do not eat, drink, smoke or apply cosmetics where immunodiagnostic materials are being handled.
- Any equipment directly in contact with specimens and reagents as well as washing solutions should be considered as contaminated products and treated as such.
- Wear lab coats and disposable gloves when handling reagents and samples and thoroughly wash your hands after handling them.
- Avoid spilling samples or solutions containing samples.
- Avoid any contact of the TMB-Substrate and the Stopping Reagent with the skin and mucosa.
- Provide adequate ventilation.
- All materials contacted with specimens or reagents, including liquid and solid wastes, should be inactivated by validated procedures (autoclaving or chemical treatment) and disposed in accordance with applicable local law regulations.



Warning!

Conjugate, Calibrators 0-5, Control Serum contain ProClin 300.

H317: May cause an allergic skin reaction.

P261: Avoid breathing vapours.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

**Instructions for use abia AFP
AB Diagnostic Systems GmbH**



Danger!

Stopping Reagent contains 0.2M sulfuric acid.

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303 + P361 + P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/ physician.

VII. MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED WITH THE KIT:

- Distilled or deionized water.
- Automatic or semiautomatic, adjustable or preset single-channel and multi-channel pipettes with a changeable volume for a set of liquids.
- Disposable pipette tips.
- Automatic microplate washer.
- Microplate reader equipped with 450 and 405-415 nm filter.
- Open type automated analyzer with 450 and 405-415 nm filter (for automated procedure).
- Laboratory clock.

VIII. COLLECTION AND HANDLING OF SPECIMENS

Blood samples should be collected according to the current practices. Serum only may be used. Separate serum as soon as possible to avoid any hemolysis. Extensive hemolysis may affect test performance. Specimens with observable particulate matter should be clarified by centrifugation prior to testing. Suspended fibrin particles or aggregates may yield falsely positive results. Do not heat the samples. For accurate comparison to established normal values, a fasting morning serum sample should be obtained.

Samples can be stored at 2-8 °C not more than for 48 hours; they should be frozen at or below at -20 °C. Avoid repeated freeze/thaw cycles. Samples that have been frozen and defrosted more than 1 time cannot be used. Samples with expressed with expressed bacterial growing, hemolysis, hyperlipidemia and which were preserved by sodium azide must not be analyzed.

IX. PREPARATION OF THE REAGENTS

1. Ready to use reagents:

- **Anti-AFP-coated microtiter wells.** Each 12-strips plate (breakable wells) is wrapped in a sealed foil-lined bag. Open the bag and remove the plate. Select the number of strips/wells required for the assay. Place the unused strips/wells back into the foil-lined bag; reseal the foil-lined bag in a Zip-Lock plastic bag. Do not remove desiccant.
- **Calibrators 0-5;**
- **Control Serum;**
- **Conjugate;**
- **TMB-Substrate;**
- **Stopping Reagent.**

2. Reagents to prepare:

- **Working Washing Solution.** Thoroughly shake Washing Solution concentrate. To make Working Washing Solution take required amount of concentrate and mix with distilled or deionized water (1:24 ratio) in a separate vial. Thoroughly mix the solution. Thoroughly mix the solution.

X. TEST PROCEDURE

Note: Before use, allow reagents to reach room temperature for 30 min.

1. Format the microplate wells for each Calibrators, Control Serum and patient specimens to be assayed in duplicate, add one or two wells for TMB control (blank).

2. Pipette 100 µl Conjugate into each well, except blank. The Conjugate should be added to the wells directly before addition of the control specimens and tested specimens.

3. Pipette 25 µl of each Calibrators, Control Serum and samples with new disposable tips into appropriate wells. And then, mix well using pipette tip. Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.

4. Swirl the microplate gently for 30 seconds after adding of samples and Conjugate to mix, cover the strips with a protective film and incubate for 120 minutes at room temperature (here 20-25 °C).

5. Aspirate the contents of the wells into the container with disinfecting solution. Wash the wells 5 times with 300µl of Working Washing Solution per well and remove Working Washing Solution using a washer into the container with disinfecting solution. Tap the plate firmly against absorbance paper to ensure that it is dry – the residual volume must be lower than 10 µl (the use of a washer is recommended). Do not allow the microplate to dry between the end of the washing operation and the reagent distribution.

6. Pipette 100 µl of TMB-Substrate into each well.

7. Incubate for 30 minutes at room temperature in a dark place.

8. Pipette 150 µl of Stopping Reagent into each well. Gently mix for 5-10 seconds.

9. Read the plate on microplate reader at 450 nm. In case of overflow absorbance values, read at 405-415 nm. Reading must be completed within 20 minutes after addition of the Stopping Reagent.

Scheme of the assay is represented in Annex.

Spectrophotometric verification of reagent pipetting

The presence of Conjugate + sample in the well can be verified by automatic reading at 405 nm. Each well containing sample and Conjugate must have an OD higher than 0.100.

10. Automated analyzers

Validated test protocols and dilution tables of reagent working solutions for different EIA-analyzers can be obtained from the manufacturer upon request (see section XIV). For the instrumentation without established validated protocol follow the section “TEST PROCEDURE” and ensure all requirements described in the section “PRECAUTIONS” are fulfilled. All protocols for automated analyzers must be fully validated before use.

When preparing working reagent solutions for automated EIA procedure, it is necessary to consider “dead” volume of vials and containers used for loading working solutions in the EIA analyzer.

XI. CALCULATION OF RESULTS

1. Calculate the mean optical density of each Calibrator duplicate.
2. Calculate the mean optical density of each unknown duplicate.
3. Subtract the mean absorbance value of the “blank” from the mean absorbance values of the Calibrators, Control Serum and serum samples.
4. Draw a calibration curve on linear graph paper with the mean optical densities on the Y-axis and the calibrator concentrations on the X-axis. If immunoassay software is being used, a 4-parameter curve is recommended.
5. In case of overflow absorbance values at 450 nm, read the results at 405-415 nm.
6. If a sample reads more than value of Calibrator 5 then dilute it with Calibrator 0. The result obtained should be multiplied by the dilution factor.

Typical tabulated data

Calibrator	OD1	OD2	Mean OD-blank	Value (IU/ml)
0	0.041	0.042	0.000	0
1	0.094	0.091	0.051	5
2	0.274	0.268	0.230	20
3	0.647	0.625	0.595	50
4	1.773	1.733	1.712	150
5	2.613	2.675	2.603	300
Unknown	0.515	0.510	0.471	40

This data is for the purpose of illustration only, and should not be used to calculate unknowns. Each user should obtain his or her own data and standard curve.

Test Validation

For the test to be valid the following criteria must be met. If these criteria are not met the test should be considered invalid and should be repeated.

1. **Blank OD:** The absorbance value should not be more than 0.1 at 450 nm.
2. The absorbance (OD) of **Calibrator 5** should not be less than 1.3 at 450 nm.
3. Calculated Value of **Control Serum** should be within established range.

XII. PERFORMANCE CHARACTERISTICS OF abia AFP

1. Assay Dynamic Range

The range of the assay is between 0-300 IU/ml.

2. Analytical sensitivity

The lower detection limit is calculated from the standard curve by determining the resulting concentration of the mean OD of Calibrator 0 (based on 24 replicate analyses) plus 2 SD.

Therefore, the sensitivity of the abia AFP kit does not exceed **1 IU/ml**.

3. Specificity (cross-reactivity)

No cross-reactivity was found for triglycerides up to 661 mg/dl, for cholesterol up to 204 mg/dl, for E.coli Ab, for rheumatoid factor up to 77.4 U/ml, for hemoglobin up to 2650 mg/l, for bilirubin up to 502 µmol/l, for HAMA up to 507.1 ng/ml.

4. Precision

Intra-assay precision

The within assay variability is shown below:

Sample	n	Mean, ng/ml	SD	CV, %
1	27	145.04	6.19	4.27
2	27	41.18	1.08	2.63

Inter-Assay precision

The between assay variability is shown below:

Sample	n	Mean, ng/ml	SD	CV, %
1	3	145.08	9.39	6.47
2	3	40.08	2.43	6.07

5. Recovery

Samples were prepared by adding defined amounts of AFP (Calibrator 2) to Control Serum. The results (in IU/ml) are tabulated below:

Sample	Number of repeats	Measured concentration, IU/ml	Expected concentration, IU/ml	Recovery, %
Serum sample	9	41.08	-	-
Calibrator 2	3	20.00	-	-
Serum sample + Calibrator 2	9	32.21	30.54	105.48%

6. Linearity

The high-AFP concentration serum sample (H) was diluted with Calibrator 0 which was used as the low-concentration sample (L). The results (in IU/ml) are tabulated below:

Dilution	Coded Concentration (IU/ml)	Replicates	Mean	Recovered, %
L	-	3	0.00	-
0.975L+0.025H	8.41	3	8.29	98.48
0.95L+0.05H	16.83	3	16.03	95.28
0.9L+0.1H	33.66	3	34.40	102.21
0.8L+0.2H	67.32	3	63.36	94.12
0.7L+0.3H	100.97	3	107.79	106.75
0.6L+0.4H	134.63	3	142.20	105.62
0.5L+0.5H	168.29	3	160.94	95.63
0.4L+0.6H	201.95	3	190.19	94.18
0.3L+0.7H	235.60	3	230.55	97.85
0.2L+0.8H	269.26	3	266.42	98.94
0.1L+0.9H	302.92	3	298.06	98.40
H	-	3	336.58	

7. Expected normal Value

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

Group	Number of samples	Range, IU/ml	Median value, IU/ml	Claimed reference limit
Adult	131	0-6	2.2	<10

Unit Conversion Calculator: IU/ml x 1.205=ng/ml; ng/ml x 0.83=IU/ml

8. High dose hook effect

The assay was tested for a high dose hook effect. Up to a AFP concentration of 1500 IU/ml no hook effect was observed.

9. Accuracy

The abia AFP kit was compared with a Chemiluminescent microparticle immunoassay as a reference test. The total number of specimens was 40. The values ranged from 1.4 to 197 IU/ml. The least square regression equation and the correlation coefficient were computed for abia AFP in comparison with the reference method. The least square regression analysis was $y=0.9459(x) + 2.822$ with correlation coefficient 0.97.

XIII. LIMITS OF THE TEST

- All the reagents within the kit are calibrated for the direct determination of AFP in human serum. The kit is not calibrated for the determination of AFP in saliva, plasma or other specimens of human or animal origin.
- Any improper handling of samples or modification of this test might influence the results.
- Only Calibrator 0 may be used to dilute any high serum samples. The use of any other reagent may lead to false results.
- The results obtained with this kit should never be used as the sole basis for clinical diagnosis. Any laboratory result is only a part of the total clinical picture of the patient.
- Some individuals may have heterophilic antibodies to mouse or other animal proteins that can possibly interfere in this assay. Therefore, the results from any patients who have received preparation of mouse antibodies for diagnosis or therapy should be interpreted with caution.
- Not intended for newborn screening.

XIV. CONDITIONS OF STORAGE AND TRANSPORTATION

- **Expiry date is indicated on the packaging.** Storage and transportation conditions for the kit, conditions and terms of storage for working solutions and unused reagents are specified in table 2.
- Transportation should be done by covered transport at specified temperature in accordance with established transportation regulations. Kits transported at improper temperature cannot be used.
- Kits stored improperly cannot be used.

1	Storage conditions		
	Keep in a dark dry place at 2-8 °C. Freezing is prohibited.		
2	Transportation conditions		
	at 2-8 °C		
	at 9-25 °C	not more than during ten (10) days	
3	Conditions and terms of storage for working solutions		
	Keep in a dark dry place and in a chemically neutral vial.		
	Working Washing Solution	at 2-8 °C at 18-24 °C	For up to 28 days For up to 14 days
4	Conditions and terms of storage of unused reagents after opening		
	Keep in a dark dry place at 2-8 °C.		
	Anti-AFP-coated microtiter wells	Place the unused strips/wells back into the bag, reseal the foil-lined package in Zip-Lock plastic bag. Do not remove desiccant.	Until the kit expiration date.
	Calibrators 0-5, Control Serum, Washing Solution, Stopping Reagent	Close the vials tightly with screw caps and stored them in the manufacturer's package.	Until the kit expiration date.
Conjugate, TMB-Substrate	Close the vials tightly with screw caps and stored them in the manufacturer's package.	For two months.	



AB Diagnostic Systems GmbH

Sportfliegerstraße 4, Berlin, 12487, Germany

















Tel. +49 30 208987160, Fax: +49 30 208987199

E-Mail: info@ab-ds.de, www.ab-ds.de

XV. REFERENCES

1. Bates S.E. Clinical applications of serum tumor markers. *Ann Intern Med* 1991; 115:623-8.
2. Wu J.C., Lee S.D., Hsaio K.J., et al. Mass screening of primary hepatocellular carcinoma by alpha-fetoprotein in a rural area of Taiwan-a dried blood spot method. *Liver* 1988; 8:100-4.
3. Lee H.-S., Chung Y.H., Kim C.Y. Specificities of serum alpha-fetoprotein in HBsAg+ and HBsAg- patients in the diagnosis of hepatocellular carcinoma. *Hepatology* 1991; 14:68-72.
4. Di Bisceglie A.M., Rustgi V.K., Hoofnagle J.H., Dusheiko G.M., Lotze M.T. Hepatocellular carcinoma. *Ann Intern Med* 1988; 108:390-401.
5. Sato Y., Nakata K., Kato Y., et al. Early recognition of hepatocellular carcinoma based on altered profiles of alpha-fetoprotein. *N Engl J Med* 1993; 328:1802-6.
6. Deutsch H.F. Chemistry and biology of alpha-fetoprotein. *Adv Cancer Res* 1991; 56:253-312.
7. Inczedy J., Lengyel T., Ure A. Compendium of Analytical Nomenclature (definitive rules 1997) - The Orange Book, 3rd edition, International Union of Pure and Applied Chemistry Blackwell Science. 1998, 0-632-05127-2.

XVI. EXPLANATION OF SYMBOLS

	CE marking (European directive 98/79/CE on in vitro diagnostic medical devices)		Storage temperature limitation
	Manufacturer		Consult Instruction for use
	Date of manufacture CCYY-MM		For in vitro diagnostic use
	Expiry date CCYY-MM-DD		Sufficient for
	Batch code		Symbol “exclamation mark”
	Catalog number	Warning!	Signal word
	Fragile, handle with care		Symbol “corrosion”
	Keep away from sunlight	Danger!	Signal word
	Keep dry		Top

Scheme of the assay

1	Add	100 µl of Conjugate into each well, except blank
2	Add	25 µl of Calibrators, Control Serum in duplicates; 25 µl of samples in duplicates; one or two wells for TMB control (blank)
3	Mix	30 seconds
4	Incubate	120 min, at 20-25 °C
5	Wash the plate	Working Washing Solution, 300 µl, 5 times
6	Add	100 µl of TMB-Substrate to all wells
7	Incubate	30 min, at room temperature in a dark place
8	Add	150 µl of Stopping Reagent to all the wells
9	Mix	5-10 seconds
10	Read the optical density	450 nm, 405-415 nm

I IU/ml x 1.205=ng/ml

ng/ml x 0.83=IU/ml